# Pressure Support Ventilation Advisory System Provides Valid Recommendations for Setting Ventilator

Steven Bonett RRT, Michael J Banner PhD RRT, Neil R Euliano PhD, Carl W Peters MD, A Joseph Layon MD, and Andrea Gabrielli MD

BACKGROUND: Pressure support ventilation (PSV) should be applied so that the inspiratory muscles are unloaded appropriately. We developed a computerized advisory system that assesses the load on the inspiratory muscles to spontaneously inhale, reflected by the automatically and noninvasively measured work of breathing per minute, and tolerance for that load, reflected by spontaneous breathing frequency and tidal volume, in a fuzzy-logic algorithm that provides recommendations for setting PSV. We call this a load and tolerance strategy for determining PSV. METHODS: In a clinical validation study, we compared the recommendations from our PSV advisory system to the recommendations of experienced critical-care Registered Respiratory Therapists (RRTs) for setting PSV in patients with respiratory failure. With 76 adult patients in a university medical center surgical intensive care unit receiving PSV, a combined pressure/flow sensor, positioned between the endotracheal tube and patient Y-piece, sent measurements to the PSV advisory system. We compared the advisory system's recommendations (increase, maintain, or decrease the pressure support) to the RRTs' recommendations at the bedside, RESULTS: There were no significant differences between the RRTs' and the advisory system's recommendations (n = 109) to increase, maintain, or decrease PSV. The RRTs agreed with 91% of the advisory system's recommendations (kappa statistic 0.85, P < .001). The advisory system was very good at predicting the RRTs' pressure support recommendations ( $r^2 = 0.87, P < .02$ ). CONCLUSIONS: A load and tolerance strategy with a computerized PSV advisory system provided valid recommendations for setting PSV to unload the inspiratory muscles, and the recommendations were essentially the same as the recommendations from experienced critical-care RRTs. The PSV advisory system operates continuously and automatically and may be useful in clinical environments where experts are not always available. Key words: work of breathing; mechanical ventilation; pressure support ventilation; acute respiratory failure; respiratory monitoring; artificial neural networks. [Respir Care 2011;56(3):271–277. © 2011 Daedalus Enterprises]

### Introduction

A goal of pressure support ventilation (PSV) is to appropriately unload the inspiratory muscles of patients with compromised pulmonary function and respiratory failure. The inspiratory muscles become loaded due to decreased respiratory system compliance and/or increased respiratory system resistance. The amount of inspiratory positive pressure assist applied with PSV is typically determined by observing the patient's spontaneous breathing pattern. The inspiratory muscle load, which can be assessed by work of breathing per minute (WOB/min), is thought to be appropriate when the spontaneous breathing frequency (f) and tidal volume ( $V_T$ ) are in clinically acceptable ranges,

a comfortable breathing pattern is present, and accessory inspiratory muscles (eg, sternocleidomastoid) are not active. Because inspiratory muscle load and breathing-pattern data do not always correlate when applying PSV,<sup>2,3</sup>

SEE THE RELATED EDITORIAL ON PAGE 359

another approach is to consider using measured values for WOB/min, f, and  $V_T$ , and then using a portable computerized bedside monitor to assess the interrelationships of these values and formulate recommendations for setting PSV.

Changes in the patient's pulmonary mechanics and peak inspiratory flow rate may require a change in PSV, but highly trained, skilled clinicians may not always be available to make the assessments to properly set PSV. For these situations we propose an open-loop, computerized PSV advisory system that automatically assesses the breathing load on the inspiratory muscles (WOB/min) and tolerance for that load (f and V<sub>T</sub>) and provides recommendations for increasing, maintaining, or decreasing PSV.

Previously, we performed a multi-site clinical validation study in which we compared the recommendations of critical-care attending physicians to the recommendations of our PSV advisory system for setting PSV in patients with respiratory failure.4 The PSV advisory system provided valid recommendations for setting PSV to appropriately unload the inspiratory muscles, which were as good as the clinical judgment of attending physicians. Another group of clinicians who may make recommendations for setting PSV is critical-care Registered Respiratory Therapists (RRTs). In some hospitals RRTs are delegated, by authority of an attending physician, to make recommendations and changes to PSV settings within a specific range, and to wean patients from ventilatory support, as articulated in RRT-driven weaning protocols, for example. In other situations, experienced critical-care RRTs may be tasked to guide ventilatory management once a management algorithm is agreed upon between RRT leadership and the intensive care unit (ICU) attending medical staff. Thus, we conducted this follow-up study to compare the clinical judgment of critical-care RRTs in setting PSV to the recommendations from the PSV advisory system. We hypoth-

Steven Bonett RRT, Michael J Banner PhD RRT, Carl W Peters MD, A Joseph Layon MD, and Andrea Gabrielli MD are affiliated with the Department of Anesthesiology, University of Florida College of Medicine, Gainesville, Florida. Steven Bonett RRT is also affiliated with Cardiopulmonary Services, Shands Hospital, Gainesville, Florida. Carl W Peters MD, A Joseph Layon MD, and Andrea Gabrielli MD are also affiliated with the Department of Surgery, University of Florida College of Medicine, Gainesville, Florida. A Joseph Layon MD is also affiliated with the Department of Medicine, University of Florida College of Medicine, Gainesville, Florida. Neil R Euliano PhD is affiliated with Convergent Engineering, Gainesville, Florida.

The authors have disclosed a relationship with Respironics and Convergent Engineering.

Steven Bonett RRT presented a version of this paper at the 54th International Respiratory Congress of the American Association for Respiratory Care, held December 13-16, 2008, in Anaheim, California.

Correspondence: Michael J Banner PhD RRT, Department of Anesthesiology, University of Florida College of Medicine, Box 100254, Gainesville FL 32610. E-mail: mbanner@anest.ufl.edu.

DOI: 10.4187/respcare.00656

Table 1. Patient Data\* at Enrollment (n = 76)

Age (y)	$57.5 \pm 16$
Male/female (no.)	47/29
Weight (kg)	$88 \pm 31$
Pressure support (cm H <sub>2</sub> O)	$12 \pm 4$
PEEP (cm H <sub>2</sub> O)	$7 \pm 4$
IMV rate (breaths/min)	$2 \pm 1$
$F_{IO_2}$	$0.37 \pm 0.05$
f (breaths/min)	$19 \pm 7$
$V_{T}$ (mL/kg)	$7.0 \pm 1.2$
$\dot{V}_{E}$ (L/min)	$13.2 \pm 2.7$
P <sub>ETCO2</sub> (mm Hg)	$38 \pm 9$
$S_{pO_2}$ (%)	$97 \pm 2$
C <sub>RS</sub> (L/cm H <sub>2</sub> O)	$0.06 \pm 0.04$
$R_{RS}$ (cm $H_2O/L/s$ )	$8.5 \pm 3.8$

<sup>\*</sup> Values are mean ± SD except for male/female values.

IMV = intermittent mandatory ventilation

f = spontaneous breathing frequency

V<sub>T</sub> = spontaneous tidal volume

 $\dot{V}_E$  = spontaneous exhaled minute ventilation

 $P_{ETCO_2}$  = end-tidal partial pressure of carbon dioxide

 $C_{RS}$  = respiratory system compliance

R<sub>RS</sub> = respiratory system resistance

esized that the PSV advisory system would supply PSV recommendations comparable to those chosen by the RRTs.

### Methods

Approval was obtained from the institutional review board at the University of Florida at Shands Hospital to conduct the study in the surgical ICU. We enrolled 76 adult patients intubated or tracheostomized (tube inner diameter range 6.0-8.5 mm) with respiratory failure from various etiologies (pneumonia, pulmonary edema, sepsis, congestive heart failure, and subarachnoid hemorrhage) in the maintenance and weaning phases of ventilatory care (Table 1). Nine patients had direct lung injury from penetrating blunt chest trauma, and 6 had COPD. All the patients were breathing spontaneously, receiving PSV, receiving sedation and analgesia as needed, and hemodynamically stable. Excluded were patients who were hemodynamically unstable (eg, mean arterial blood pressure range 40–80 mm Hg, with irregular changes in heart rate and/or arrhythmias); heavily sedated (eg, high level of opioids, causing substantial changes in f and V<sub>T</sub>); or with irregular breathing patterns (eg, closed-head injury with sudden increases and decreases in f and  $V_T$ ).

All the patients were ventilated with the same type of ventilator (840, Puritan-Bennett, Pleasanton, California), set in the PSV and PEEP modes, which we combined with synchronized intermittent mandatory ventilation (SIMV) set at 2 breaths/min and a  $V_T$  of 6–8 mL/kg ideal body weight. It is a patient-care policy in our surgical ICU to

employ SIMV at a minimum of 2 breaths/min when using PSV and PEEP, and our institutional review board required that we conform to all existing patient-care policies, so we used SIMV, with a decelerating inspiratory flow waveform. For PSV we set the "% rise" setting in the range 60–80%. This setting controls the ventilator's flow rate output during inhalation and affects the rate of pressure rise (ie, the interaction of the ventilator's flow rate output and patient's inspiratory flow rate demand affect the rate of pressure rise or shape of the pressure profile during a pressure-supported breath). With all patients we set the expiratory sensitivity setting (breath-termination criterion) at 25%.

Two components constitute the open-loop PSV advisory system: a commercially available respiratory monitor (NICO, Respironics, Wallingford, Connecticut), and a laptop computer. A combined pressure/flow/carbon dioxide sensor, positioned between the endotracheal tube and the patient Y-piece, directs data to the respiratory monitor for measurement of f,  $V_T$ , exhaled minute ventilation, PSV, PEEP, and end-tidal partial pressure of carbon dioxide ( $P_{\rm ETCO_2}$ ) (Fig. 1). Compliance and resistance of the respiratory system were calculated with the least-squares method.  $^6$  S $_{\rm PO_2}$  readings were collected with a finger sensor.  $P_{\rm ETCO_2}$  and S $_{\rm PO_2}$  data were monitored to ensure that the patients were appropriately ventilated and oxygenated.

Data from the monitor go to the laptop computer where the PSV advisory system software (Convergent Engineering, Gainesville, Florida) operates an artificial neural network for immediate calculation of noninvasively measured WOB/min (WOB<sub>N</sub>/min)<sup>3</sup> and a fuzzy-logic inference system that processes the WOB<sub>N</sub>/min, f, and V<sub>T</sub> information to assess the load on the inspiratory muscles to inhale  $(WOB_N/min)$  and the tolerance for that load (f and  $V_T$ ) and provides recommendations to increase, maintain, or decrease the pressure support. We refer to this as a load and tolerance strategy (see Fig. 1). Fuzzy logic is a process of using probability distributions instead of simple "yes/no" decisions, as in a simple rule-based system, to drive the ventilator decision making.<sup>7,8</sup> The WOB<sub>N</sub>/min data were highly correlated (r = 0.91, P < .02) with invasively measured WOB/min (from esophageal pressure measurements) and were considered a very good predictor of invasively measured WOB/min in adults treated with PSV.<sup>3</sup> WOB<sub>N</sub>/min reflects the total load on the inspiratory muscles, which includes the elastic load of the respiratory system and resistive loads of the airways, endotracheal tube, and ventilator apparatus. A limitation of WOB<sub>N</sub>/min is that it cannot differentiate these component loads.

Artificial neural networks—a branch of artificial intelligence—have been used for medical applications<sup>9-12</sup> (eg, breathing-pattern recognition during spontaneous breathing and PSV).<sup>11</sup> Artificial neural networks are a class of mathematical models that are biologically inspired com-




Fig. 1. Open-loop pressure support ventilation (PSV) advisory system. Data from the pressure and flow sensors at the Y-piece go to the respiratory monitor, then to the PSV advisory system software, which also collects data from the ventilator. The PSV advisory system's artificial neural network software calculates the noninvasively measured work of breathing per minute (WOB<sub>N</sub>/min), and the fuzzy-logic inference system software uses the WOB<sub>N</sub>/min, spontaneous breathing frequency (f), and tidal volume (V<sub>T</sub>) to form its recommendation for setting PSV to appropriately unload the inspiratory muscles. In this load and tolerance strategy the load parameter, WOB<sub>N</sub>/min, reflects the load on the inspiratory muscles during spontaneous inhalation, and the tolerance parameters, f and V<sub>T.</sub> reflect the tolerance for the load on the inspiratory muscles. In this figure the PSV advisory system determined that the current pressure support setting of 10 cm H<sub>2</sub>O is insufficient because WOB<sub>N</sub>/min is too high (17 J/min) and tolerance for that load is borderline (f is 32 breaths/min and V<sub>T</sub> is low at 5 mL/kg ideal body weight), so the system recommends increasing PSV.

puter programs designed to simulate the way in which the human brain processes information. 13,14 An artificial neural network can detect complex nonlinear relationships between dependent and independent variables in data that a human brain may fail to detect. They acquire knowledge by detecting patterns and relationships to facilitate learning; they are trained through experience. Basically, an artificial neural network learns from data and captures the knowledge contained in the data. An artificial neural network is formed by hundreds of single units (artificial neurons, or processing elements) connected with coefficients or weights (analogous to synaptic connections in the cen-

tral nervous system), which constitute the neural system, and are organized into layers. The power/intelligence of neural computations is achieved from connecting neurons in a network. During training, the inter-unit connections are optimized until the errors in predictions are minimized and the artificial neural network reaches a specified level of accuracy. Once the network is trained, it can be given new input information to predict the output. A "multi-layer perceptron" (the most common type of artificial neural network) can be thought of as a combination of the weighted sum of a number of logistic regression models. Each logistic regression model (a "perceptron" in neural network terminology) has its own set of weights, and its output is weighted and added to the outputs of the other perceptrons to create the final prediction. WOB/min is the predicted variable of the artificial neural network in our PSV advisory system.

The training of our artificial neural network was previously described.<sup>3</sup> Five input elements (predictors) are used:

- Spontaneous minute ventilation (not including SIMV breaths) correlates directly with WOB/min.
- Increased intrinsic PEEP is associated with increased WOB/min, and vice versa.
- Lower inspiratory trigger pressure depth just before the ventilator triggers during PSV is associated with increased WOB/min, and vice versa.
- Lower inspiratory flow rise time during a PSV breath is associated with increased WOB/min, and vice versa.
- Higher respiratory muscle pressure (pressure = (V<sub>T</sub>/respiratory system compliance) + (inspiratory flow rate × respiratory system resistance)) is associated with increased WOB/min, and vice versa.<sup>15,16</sup>

One goal of the PSV advisory system is to maintain  $WOB_N/min$ , f, and  $V_T$  in the green ranges shown in Figure 1, which were derived, in part, from normal WOB values,  $^{3,15}$  a WOB outcome study,  $^2$  and our clinical experience.  $^3$  In adults the normal range of WOB/min is 4-8 J/min.  $^{15}$  Based on 15 years of experience measuring WOB in patients receiving ventilatory support, adults treated with PSV tolerate a maximum WOB/min of up to about 10-12 J/min. WOB/min values greater than 12-15 J/min in adults are not well tolerated and require a higher pressure support level. In adults treated with PSV it was found that WOB of 2-8 J/min was well tolerated.  $^3$  Another clinical study reported that WOB maintained in a fairly normal range during PSV was well tolerated.

Our PSV advisory system's target f and  $V_T$  ranges are congruent with those applied to spontaneously breathing adults receiving PSV: 10-25 breaths/min, and 6-8 mL/kg, respectively.<sup>1,17</sup> It is reasonable to apply PSV so that inappropriately low f (eg, 4-6 breaths/min) or high f (eg, 30-40 breaths/min) and low  $V_T$  (eg, <4 mL/kg) or high

 $V_T$  (eg, > 12 mL/kg) do not occur. The advisory system treats the wider WOB<sub>N</sub>/min, f, and  $V_T$  ranges (the yellow areas in Figure 1) as cautionary. The advisory system assesses the probability of inspiratory muscle load (WOB<sub>N</sub>/min) and tolerance (f and  $V_T$ ) parameters varying in relation to one another in all ranges, and then provides an appropriate recommendation to increase, maintain, or decrease PSV. The system uses a 5-minute averaging filter for updating its recommendation.

PSV was set initially as ordered by the attending physicians. Subsequently, the RRTs' pressure support recommendations were permitted during the study. At periodic intervals the RRTs were asked to evaluate the patient to determine if PSV should be increased, maintained, or decreased. Their recommendations were based on using the traditional approach of assessing spontaneous breathing-pattern data (ie, f 15–25 breaths/min, V<sub>T</sub> 6–8 mL/kg ideal body weight), no sternocleidomastoid muscle contraction, and appearance of comfortable breathing. The RRTs were blinded to the measured respiratory data and the advisory system's recommendations.

PEEP and  $F_{\rm IO_2}$  settings, determined by the attending physicians (see Table 1), were held constant during the study period, which was over an 8-hour work shift. There were no clinically important differences in PEEP,  $P_{\rm ETCO_2}$ , or  $S_{\rm pO_2}$  at baseline or during the study, and no clinically important variations in hemodynamic parameters in any patients throughout the study.

We assessed the relationship between PSVs set by the clinicians and PSVs recommended by the advisory system. When the RRTs recommended an increase or decrease in PSV, generally the change was by 5 cm H<sub>2</sub>O from the current setting. When the advisory system recommended "Increase PSV" or "Decrease PSV," this was interpreted to mean a change of 5 cm H<sub>2</sub>O up or down, respectively, from the current setting. With this understanding, the RRTs' recommended PSV settings were regressed with the recommended settings from the PSV advisory system.

We analyzed the data with chi-square, regression analysis, and the kappa statistic.<sup>18</sup> Alpha was set at .05 for statistical significance.

# Results

PSV ranged from 5 cm  $\rm H_2O$  to 20 cm  $\rm H_2O$  during the study. The advisory system made 109 recommendations to increase, maintain, or decrease PSV. There were no clinically important differences between the RRTs' and the advisory system's recommendations (chi-square = 0.49, P = .78) (Table 2). The mean total agreement between all the RRTs' and the advisory system's recommendations was 91% (kappa statistic  $\pm$  standard error = 0.85  $\pm$  0.05, P < .001). A kappa statistic value of 0.81–1.0 is considered almost perfect agreement. There were about 3 times

Table 2. Recommendations Made by the Ventilator Advisory System, Compared to Recommendations Made by the RRTs for Setting PSV\*

RRT Recommendations	Advisory System Recommendations			
	Decrease PSV	Maintain PSV	Increase PSV	Total
Decrease pressure	65	4	0	69
Maintain pressure	6	21†	0	27
Increase pressure	0	0	13‡§	13
Total	71	25	13	109

<sup>\*</sup> Chi-square analysis. Percent agreement 91 % (99 /109), percent disagreement 9 % (10/109). A total of 109 recommendations were made by the advisory system and critical care Registered Respiratory Therapists (RRTs) to increase, maintain, or decrease pressure support ventilation (PSV). No significant differences in recommendations were made by the ventilator advisory system, compared to those of RRTs for setting PSV (chi-square = 0.49, P=.78). Overall, there was a 91% level of agreement between the recommendations of the advisory system and the RRTs. The strength of this agreement was  $K=0.85\pm0.05$  (SE) (P<.001). There were about 3 times as many recommendations to decrease PSV, compared to maintaining it († P<.01), about 5 times as many recommendations to decrease PSV than increase it ( $\ddagger P<.011$ ), and almost twice as many recommendations to maintain PSV rather than increase it (\$ P<.02).

as many recommendations to decrease PSV as to maintain it (P = .01), about 5 times as many recommendations to decrease PSV as increase it (P = .01), and about twice as many recommendations to maintain PSV as to increase it (P = .02) (see Table 2).

There was an excellent positive correlation between PSV ordered and that recommended by the advisory system:  $r=0.93,\,P<.02,\,r^2=0.87.$  The advisory system predicted or explained 87% of the variance in the RRTs' recommendations (Fig. 2).

The WOB<sub>N</sub>/min values range was 1.2–18.6 J/min. The RRTs recommended "Maintain PSV" (implying appropri-




Fig. 2. Relationship between the Registered Respiratory Therapists' (RRTs') recommendations and the automated advisory system's recommendations on adjusting the pressure support level during pressure support ventilation (PSV). The clusters of data points (eg, at 5, 10, 12, and 15 cm  $H_2O$ ) represent considerable overlap of many of the recommendations at those levels. There was an excellent and significant correlation between the RRTs' and the PSV advisory system's recommendations (r = 0.93), and the advisory system was a very good predictor of the pressure-support level recommended by the RRTs ( $r^2 = 0.87$ ).

ate inspiratory muscle loading) 27 times, and the advisory system agreed 25 times (93% agreement). Under that condition the mean  $\pm$  SD WOB<sub>N</sub>/min was 6.0  $\pm$  4 J/min.

## Discussion

An inspiratory muscle load (WOB<sub>N</sub>/min) and tolerance (f and V<sub>T</sub>) strategy in an open-loop, computerized PSV advisory system provided valid recommendations for setting PSV in adults with respiratory failure. The advisory system's recommendations were congruent with those of the critical-care RRTs evaluating patients at the bedside 91% of the time (see Table 2). The RRTs recommended decreasing PSV by a 3-to-1 margin compared to maintaining PSV, and by a 6-to-1 margin compared to increasing PSV, and the advisory system's recommendations were in accordance with that recommendation trend. Our load and tolerance strategy appears to be appropriate in view of the high level of agreement between the recommendations of the advisory system and the RRTs. At no time did the advisory system recommend a change in the opposite direction to that of the RRTs.

Concerning the 9% disagreement for all recommendations, at times there may have been philosophical differences between the thinking of some of the RRTs and the advisory system's fuzzy-logic inference system. Or some RRTs may have misinterpreted a patient's breathing pattern data and inferred a need to maintain the current PSV setting while the advisory system recommended decreasing it. Although the difference was very small, that was the most common disagreement (see Table 2). Another explanation may be a clinical decision to totally unload the inspiratory muscles with higher pressure support  $(> 20 \text{ cm H}_2\text{O})$  or to exercise them with lower pressure support (5 cm H<sub>2</sub>O). In either case, the PSV advisory system is not privy to a clinician's intentions and is incapable of "knowing" decisions to rest or exercise inspiratory muscles. This can cause disagreements, because the advisory system is unable to assimilate these types of decisions into its recommendations for setting PSV.

The load and tolerance operational strategy of our open-loop system is novel (see Fig. 1). In one respect, it determines WOB/min or power of breathing automatically, obviating esophageal balloon catheter and complex equipment.<sup>3</sup> This method greatly simplifies measurement of WOB, especially during prolonged ventilatory support, where it may be impractical to employ an esophageal balloon catheter. WOB/min reflects the load on the inspiratory muscles to spontaneously inhale, <sup>15,19</sup> while f and V<sub>T</sub> reflect a patient's tolerance for that load. Breathing load and tolerance are different but related concepts and should be used together. We believe that assessing both the inspiratory muscle load and the tolerance for that load (rather than assessing only patient tolerance) provides a more com-

prehensive understanding of a patient's ability to breathe spontaneously. Both concepts are important, and neither should be ignored. Our advisory system uses inspiratory muscle load and tolerance data in a complementary manner to determine a recommended pressure support level.

Other open-loop ventilatory advisory systems have been described for adults<sup>20,21</sup> and neonates.<sup>22</sup> With the neonatal ventilatory advisory system there was high agreement between the clinicians' and the advisory system's recommendations for setting the ventilator to treat ventilation and oxygenation abnormalities. In another study, an approach akin to an open-loop PSV system was employed, in which the response or input variable was WOB per breath, directly measured with an esophageal balloon catheter, and the output variable was PSV.<sup>2</sup> The clinicians maintained WOB per breath in a specific and fairly normal physiologic range by increasing or decreasing PSV as needed. When WOB per breath was maintained in that range with PSV, adults weaned faster than with the traditional method.

Intelligent decision-support systems for mechanical ventilation employing closed-loop control have been recommended.<sup>23</sup> Closed-loop feedback systems for controlling PSV have been described.<sup>24-26</sup> For an intelligent decision-support system closed-loop ventilator to be useful, it should be effective, safe, and easy to use at the patient's bed-side.<sup>23</sup> Based on the findings from the present study, our open-loop system meets those criteria: it was effective in recommending appropriate PSV settings, appears to be safe, in that its recommendations were validated, and because it is automatic, it is easy to use. Our open-loop system may be the basis for a closed-loop ventilator feedback system that uses WOB<sub>N</sub>/min and spontaneous f and V<sub>T</sub> data as input variables to be directed back to a ventilator and used for automatically adjusting PSV.

Our PSV open-loop advisory system differs from Smart-Care (Dräger Medical, Lübeck, Germany),25 which is a closed-loop system in that it assesses WOB/min, the load on the inspiratory muscles to spontaneously inhale, f, and V<sub>T</sub>, which reflect only the tolerance for the load. Smart-Care uses f,  $V_T$ , and  $P_{\mathrm{ETCO}_2}$  as feedback or input data to automatically regulate PSV to maintain f between 12 and 28 breaths/min,  $V_T > 0.3 L$  (or 0.25 L and if the patient weighs < 55 kg) and maintain  $P_{\rm ETCO_2} <$  55 mm Hg (or < 65 mm Hg if the patient has COPD). A potential shortcoming of that approach may be that WOB/min is not measured and taken into consideration when automatically regulating PSV. Because WOB/min is not measured, neither SmartCare nor the clinician knows whether the load on the inspiratory muscles is too high or low, so Smart-Care might apply an inappropriate amount of PSV, which could prolong weaning. This may explain, in part, why reductions in weaning duration initially demonstrated with SmartCare were not confirmed when it was compared to weaning managed by experienced clinicians.<sup>25</sup>

Our system differs also from proportional assist ventilation (PAV+, Puritan Bennett, Pleasanton, California),<sup>27</sup> which is a closed-loop system designed to automatically increase or decrease the inspiratory assist pressure in proportion to the patient's spontaneous WOB per breath, which is not measured. PAV+ uses a simple equation related to respiratory muscle pressure to calculate WOB per breath. Our approach was to employ an artificial neural network, which is a substantially more complex mathematical model that uses coefficients calculated from a clinically collected data set. A study of patients on PAV+ showed that the WOB-per-breath data displayed on the ventilator were inaccurately low, calling that method into question.<sup>28</sup> This may offer some explanation of a clinical study that reported excessive WOB when PAV+ was applied at percentage support settings lower than 30%.<sup>27</sup>

A possible limitation of the present study is that to fully validate the advisory system a similar study with a wider variety of patient conditions and performed over the total duration of ventilatory support may be needed.

## **Conclusions**

Our open-loop computerized PSV advisory system's recommendations on increasing, maintaining, or decreasing PSV were in close agreement with the RRTs' recommendations. The load and tolerance strategy in the advisory system automatically provides valid recommendations for setting PSV to appropriately unload the inspiratory muscles, that are as good as the clinical judgment of experienced critical-care RRTs. The advisory system is automatic and continuous, so in clinical situations with limited staff (eg, large ICUs, remote work areas), where skilled critical-care RRTs are not always available to perform bedside assessments, the system may be useful in recommending appropriate PSV settings. We recommend clinical outcomes studies to further evaluate this PSV advisory system, such as a study of its usefulness in patients with more difficult forms of respiratory failure.

### REFERENCES

- MacIntyre NR. Respiratory function during pressure support ventilation. Chest 1986;89(5):677-681.
- Kirton OC, DeHaven CB, Hudson-Civetta J, Morgan J, Windsor J, Civetta JM. Re-engineering ventilatory support to decrease days and improve resource utilization. Ann Surg 1996;224(3):396-404.
- Banner MJ, Euliano NR, Brennan V, Peters C, Layon AJ, Gabrielli A. Power of breathing determined non-invasively with use of an artificial neural network in patients with respiratory failure. Crit Care Med 2006;34(4):1052-1059.
- Banner MJ, Euliano NR, MacIntyre NR, Layon AJ, Bonett S, Gentile MA, et al. Ventilator advisory system employing load and tolerance

# PRESSURE SUPPORT VENTILATION ADVISORY SYSTEM

- strategy recommends appropriate pressure support ventilation settings: multi-site validation study. Chest 2008;133(3):697-703.
- MacIntyre NR, Ho L. Effects of initial flow rate and breath termination criteria on pressure support ventilation. Chest 1991;99(1): 134-137.
- Volta CA, Marangoni E, Alvisi, V, Capuzzo M, Ragazzi R, Pavanelli L, Alvisi R. Respiratory mechanics by least squares fitting in mechanically ventilated patients: application on flow-limited COPD patients. Intensive Care Med 2002;28(1):48-52.
- MacIntyre NR. Invasive mechanical ventilation in adults: conference summary. Respir Care 2002;47(4):508-518.
- Blackwood B. Commentary: Nemoto T et al (199). Automatic control of pressure support mechanical ventilation using fuzzy logic. Nurs Crit Care 2008;13(3):178-179.
- Baxt WG. Application of artificial neural networks to clinical medicine. Lancet 1995;346(8983):1135-1138.
- Rodvold DM, McLoad DG, Brandt JM, Snow PB, Murphy GP. Introduction to artificial neural networks for physicians: taking the lid off the black box. Prostate 2001;46(9):39-44.
- 11. Leon MA, Lorini FL. Ventilation mode recognition using artificial neural networks. Comp Biomed Res 1997;30(4):373-378.
- Patel JL, Goyal RK. Applications of artificial neural networks in medical science. Curr Clin Pharmacol 2007;2(3):217-226.
- Agatonovic-Kustrin S, Beresford R. Basic concepts of artificial neural network (ANN) modeling and its application in pharmaceutical research. J Pharm Biomed Anal 2000;22(5):717-727.
- Traeger M, Eberhart A, Geldner G, Morin AM, Putzke C, Wulf H, Eberhart LH. [Artificial neural networks. Theory and applications in anesthesia, intensive care and emergency medicine]. Anaesthesist 2003;52(11):1055-1061. Article in German.
- Millic-Emilli J. Work of breathing. In: Crystal RG, West JB (editors). The Lung. New York: Raven Press; 1991:1065-1075.
- Jolliet P, Tassaux D. Clinical review: patient-ventilator interaction in chronic obstructive pulmonary disease. Crit Care 2006;10(6):236-242.

- Sottiaux TM. Patient ventilator interactions during volume support ventilation: asynchrony and tidal volume instability: a report of three cases. Respir Care 2001;46(3):255-262.
- Ludbrook J. Statistical techniques for comparing measures and methods of measurement: a critical review. Clin Exp Pharm Phys 2002; 29(7):527-536.
- Otis AB, Fenn WO, Rahn HL. Mechanics of breathing in man. J Appl Physiol 1950;2(11):592-607.
- Kwok HF, Linkins DA, Mahfouf M, Mills GH. Adaptive ventilator F<sub>1O2</sub> advisor: Use of non-invasive estimates of shunt. Artif Intell Med 2004;32(3):157-159.
- Rutledge GW, Thomsen GE, Farr BR, Tovar MA, Polaschek JX, Beinlich IA, et al. The design and implementation of a ventilatormanagement advisor. Artif Intell Med 1993;5(1):67-68.
- Belal SY, Taktak AF, Nevill A, Spencer A. An intelligent ventilation and oxygenation management system in neonatal intensive care using fuzzy trend template fitting. Physiol Meas 2005;26(8):555-570.
- Tehrani FT, Roum JH. Intelligent decision support systems for mechanical ventilation. Artif Intell Med 2008;44(3):171-182.
- 24. Jouvet P, Farges C, Hatzakis G, Monir A, Lesage F, Dupic L, et al. Weaning children from mechanical ventilation with a computer driven system (closed-loop protocol): a pilot study. Pediatr Crit Care Med 2007;8(5):425-432.
- Rose L, Presneill JJ, Johnston L, Cade JF. A randomized, controlled trial of conventional versus automated weaning from mechanical ventilation using Smart Care/PS. Intensive Care Med 2008;34(10): 1788-1795.
- Pauldine R, Beck G, Salinas J, Kaczka DW. Closed-loop strategies for patient care systems. J Trauma 2008;64(4 Suppl):S289-S294.
- Ruiz-Ferrón F, Machado J, Morante A, Galindo S, Castillo A, Rucabado L. [Respiratory work and pattern with different proportional assist ventilation levels]. Med Intensiva 2009;33(6):269-275. Article in Spanish.
- Al-Rawas N, Euliano N, Gabrielli A, Layon AJ, Banner MJ. Inaccurate work of breathing displayed on ventilator during proportional assist ventilation. Anesthesiology 2009;111(9 Suppl):A272.

This article is approved for Continuing Respiratory Care Education credit. For information and to obtain your CRCE (free to AARC members) visit

RCJournal.com

