

Classification and Quantification of Patient-Ventilator Interactions: We Need Consensus!

In this issue of *RESPIRATORY CARE*, Vargas et al present the results of a randomized study comparing neurally-adjusted ventilatory assist (NAVA) versus variable pressure support (VPS) ventilation.¹ The goal was to evaluate the effect of both modes on synchrony and work of breathing. The study was carefully planned, and a fair amount of information is given for us to review and compare. At the end, no major differences were found. To the reader, this may seem as just another trial with evidence of no difference. However, there are many details, based on how the modes work and are tested, that can give perspective to their results.

A mode of ventilation is a predetermined pattern of patient-ventilator interaction that can be classified by 3 characteristics: (1) control variable, that is, pressure or volume; (2) breath sequence, that is, continuous mandatory ventilation, intermittent mandatory ventilation, or continuous spontaneous ventilation; and (3) targeting scheme(s).² Such taxonomy is important for distinguishing arbitrary brand names from generic classifications, just as is done in pharmacology. Once taxonomy aspects are standardized, it is easier to evaluate performance in terms of clinical outcomes (eg, mortality), technological features serving the goals of ventilation (ie, safety, comfort, and liberation),³ patient-ventilator interaction (ie, synchrony issues),⁴ and/or discrete physiological outcomes (eg, work of breathing). Vargas et al¹ focused on compared physiological outcomes (synchrony and work of breathing) between NAVA and VPS. Their study allows us to discuss a framework to compare modes. Evidently, modes differ substantially with respect to technical capabilities as well as settings used to manage them. Thus, a systematic comparison of mode characteristics could help understand differences and potential sources of bias when assessing performance. Table 1 summarizes a construct to assess the differences between NAVA and VPS with pressure support (PS) added as a comparator.

Dr Mireles-Cabodevila is a co-owner of a patent for mid-frequency ventilation. He discloses relationships with the American College of Physicians, Elsevier, and Jones & Bartlett publishers. Dr Gama de Abreu was granted a patent on variable pressure support (VPS) ventilation that is licensed to Dräger Medical. He discloses relationships with Ambu, ZOLL, Lungpacer, GE Healthcare, Dräger Medical, and Novalung.

DOI: 10.4187/respcare.10114

The first insight is that even though these modes control inspiratory pressure and have the same breath sequence, the targeting scheme (ie, feedback control software) is different. NAVA uses servo targeting and VPS biovariable targeting. NAVA uses the diaphragm electrical activity

SEE THE ORIGINAL STUDY ON PAGE 503

(EA_{di}) to control the timing and magnitude of the pressure delivered. The pressure delivered in NAVA is proportional to the amplitude of the EA_{di} signal, and the proportionality is regulated by the set NAVA level (cm H₂O/mV of EA_{di}). VPS is similar to PS in the way we set trigger, cycle and inspiratory pressure, and the average or mean PS. However, in VPS, the operator sets also the variability of the inspiratory pressure as a percentage (0–100%) of the target inspiratory pressure. The delivered inspiratory pressure values are defined randomly, but follow a gaussian distribution (extreme values are rare, whereas values around the mean occur more frequently). For example, a setting of 100% means that the inspiratory pressure for any given breath varies randomly from 100% below to 100% above (double) the set inspiratory pressure target. If the maximum support achievable is limited by the maximum pressure setting, the variability is reduced in order to not exceed that limit, while keeping a gaussian distribution of PS levels. Currently, NAVA is available only in the Getinge Servo ventilator, whereas VPS is available only in the Dräger ventilator V500. This is important because the programming and performance may change if other ventilator platforms develop or if there are updates to the ventilator software.

Table 1 helps highlight potential differences in how the ventilators may perform when evaluated. For instance, when evaluating trigger performance, NAVA uses the EA_{di} signal that leads (ie, occurs before) the airway pressure and flow signals (Fig. 1) that trigger and cycle inspiration for PS and VPS. The method to set the trigger and cycle thresholds for each mode may lead to systematic differences that can be ascribed to the operator rather than the mode per se. In the study of Vargas et al,¹ an arbitrary fixed trigger threshold (2 L/min) was used for VPS rather than the maximum sensitivity that was not associated with false triggering.⁵ Although there was no difference noted in this study, in theory, this may bias measurements of trigger delay and

Table 1. Technical Comparison of Modes

	NAVA	VPS	PS	Comparison Comments
TAG	PC-CSVr	PC-CSVb	PC-CSVs	The TAG indicates that these are similar modes that differ in their targeting schemes.
Trigger variable	<u>Adjustable</u> (EA _{di} threshold)	<u>Adjustable</u> (flow/pressure threshold)	<u>Adjustable</u> (flow/pressure threshold)	EA _{di} signal leads airway pressure/flow signal. Trigger sensitivity affects performance of mode.
Inspiratory pressure	<u>Variable</u> (proportional to EA _{di})	<u>Variable</u> (set variability around set mean inspiratory pressure)	<u>Fixed</u> (set inspiratory pressure)	Operator chosen settings (NAVA level, % variation, and set inspiratory pressure) will affect inspiratory pressure and thus performance of the mode. See targeting scheme.
Cycle variable	<u>Fixed</u> (at 70% of peak EA _{di})	<u>Adjustable</u> (% peak inspiratory flow)	<u>Adjustable</u> (% peak inspiratory flow)	EA _{di} signal leads airway pressure/flow signal. Cycle sensitivity may affect performance of mode. Neuromechanical uncoupling is not detected by airway or P _{es} signals.
Targeting scheme	<u>Servo</u> Inspiratory pressure proportional to patient effort as measured by EA _{di} .	<u>Biovariable</u> Inspiratory pressure varies breath to breath based on an equation to generate variability.	<u>Set Point</u> Inspiratory pressure is constant breath-to-breath.	Targeting schemes vary according to how well they serve the goals of ventilation (ie, safety, comfort, liberation). The comparison needs to account for these differences.

NAVA = neurally adjusted ventilatory assist
VPS = variable pressure support
PS = pressure support
TAG = taxonomic attribute grouping
PC = pressure control
CSV = continuous spontaneous ventilation
r = respiratory system resistance
b = biovariable targeting scheme
s = set-point targeting scheme¹
EA_{di} = electrical signal of the diaphragm
P_{es} = esophageal pressure

trigger work against VPS. Although this potential bias might be negligible from a clinical standpoint, differences may yield conclusions with clinical relevance, such as choosing one mode over the other.

With regard to inspiratory pressure, NAVA and VPS have different targeting schemes. In the study of Vargas et al,¹¹ the NAVA level and the set inspiratory pressure in VPS were adjusted to achieve a tidal volume of 6–8 mL kg. However, for VPS, the inspiratory pressure variability was set to 100%, which means that the patient will get a large range of inspiratory pressures around the set mean inspiratory pressure. From a synchrony standpoint, the larger variability range may predispose VPS to both underassistance and overassistance⁶ of some breaths. Overassistance might explain the significantly higher incidence of failed triggers in the study of Vargas et al.¹ On the other hand, the large range of pressure variation may also have a protective effect against continuous overassistance.

The cycling criterion is also relevant. NAVA uses a fixed threshold (70% of the peak EA_{di}), whereas in VPS the threshold can be adjusted. Vargas et al chose to use a fixed

threshold of 25% of the peak flow.¹ Although these fixed settings may lead to fair comparisons, it can also explain the nonsignificant but higher incidence of delayed cycle for VPS. In terms of NAVA, the fixed threshold is a known cause of double trigger⁶ and has been frequently reported.^{7–11}

This seems to be a technical artifact due to the fixed cycle threshold in the setting of irregular EA_{di} activation. The waveform meets the definition criteria used in most trials; that is, a double trigger is 2 breaths (2 or more triggers leading to pressurization) separated by an expiratory time less than half the mean inspiratory time.^{4,12} However, not all double-trigger events are the same. Consider the double trigger that occurs in volume control versus pressure control. In volume control, the patient will receive double the tidal volume (colloquially referred as breath-stacking), which can lead to volutrauma. In pressure control, the patient will receive the second breath at the set pressure for the set inspiratory time, and the tidal volume will depend on the respiratory system compliance and resistance as well as the pressure generated by the ventilatory muscles. The risk of volutrauma is lower; however, the




Fig. 1. Neuromechanical uncoupling and its effects on patient ventilator synchrony. EA_{di} = diaphragm electrical activity. P_{es} = esophageal pressure. P_{aw} = airway pressure.

All breaths are pressure control. Green solid line is start of neural signal; green dashed line: start of inspiratory flow from ventilator (trigger event); red solid line: peak of the EA_{di} (end of neural time); red dashed line: end of inspiratory flow from ventilator (cycle event); blue dashed lines: difference between maximum EA_{di} and minimum P_{es}. Shaded area: period of neural inspiratory time. A: Patient with overassistance. EA_{di} has low amplitude and short duration. Note delay in trigger event compared to start of EA_{di} signal. Note also the delay in cycle event well after peak EA_{di}. Airway pressure is out of synchrony with EA_{di}. B: Patient with underassistance. There is a similar trigger delay as in A, but the cycle delay is much less. However, in this case delayed cycling is evident by rise in end-inspiratory pressure as the ventilatory muscles relax.²⁰ Note the dissociation of the maximum EA_{di} and minimum P_{es}, consistent with neuromechanical uncoupling. Modified from Reference 8 and Reference 21.

theoretical risk for diaphragm injury via eccentric contractions increases. In PS, double trigger will lead to repressurization of the airway that will cycle off again per usual criteria, the effects of which are unknown but likely involve less risk of volutrauma and diaphragm injury. In NAVA, double trigger leads to brief pressurizations and volume change, the effects of which are unknown, but likely similar to PS. These differences are not accounted by definitions or by lumped indexes (more on this below).

The definitions and methods used to evaluate mode performance (in terms of esophageal pressure [P_{es}], EA_{di}, or waveform analysis) may lead us to favor one mode over another.¹³ For example, the way we measure and define synchrony. Both the EA_{di} and the P_{es} have been heralded as the accepted standard for evaluation of patient-ventilator interactions. Both are invasive and provide relevant information regarding timing and magnitude of the effort. However, the presence of neuromechanical uncoupling may affect the interpretation of the interactions when using

P_{es} versus EA_{di}. Neuromechanical uncoupling refers to the difference between neural signal and the mechanical response (change in airway flow/pressure) of the respiratory system. In Figure 1, one can note the discrepancies that would occur if the P_{es} or the EA_{di} was used to measure the discordances. Importantly, the P_{es} signals will not detect uncoupling events, which can be relevant in how we classify interactions.^{8,14} Neuromechanical uncoupling has several clinically important causes such as air flow obstruction and overassistance. The science of patient-ventilator interaction is evolving, and we need to define a set of definitions and standards for measurement and reporting of synchrony that accounts for these nuances.

Finally, we often compare modes using measures of global synchrony, such as the asynchrony index. Such an index may be calculated in different ways; some use the EA_{di},^{15,16} others the P_{es},¹⁷ and others waveform analysis.¹⁸ In all these cases, the same threshold of 10% is used, which likely has different performance across methods. The type of patient-ventilator discordance included in the global index will also affect comparisons (eg, including discordances that deal with work/demand vs those that involve timing). By lumping all discordances together, we lose precision in comparing modes. Some synchrony issues are clinically irrelevant and may simply be an analytical distraction. Some may be due to inappropriate settings. Other relevant measures, such as overassistance or early trigger (also known as reverse trigger), are not included. Any asynchrony index needs to be taken with these caveats in mind. In the study of Vargas et al,¹ modes did not differ significantly with respect to the asynchrony index, yet we can see that many of the interactions that influenced the incidence of the global index were perhaps due to the settings and to technical issues other than mode performance per se.

Patient-ventilator synchrony is a commonly evaluated outcome in mechanical ventilation. It is a flourishing area of research and clinical application. Based on recent studies, it might also affect clinical outcomes.¹⁹ It follows that our assessment of synchrony needs to be reliable and clear. At this time, we lack a unified reporting taxonomy for patient-ventilator interactions; we lack clarity on which tools to use to assess patient-ventilator interactions, and finally we lack of a standardized method to apply measuring tools. As we work to advance our understanding of mechanical ventilation, synchrony, and work of breathing, a unified multisociety consensus would be welcome. This would allow thoughtful studies such as this by Vargas et al¹ to be even more useful.

Eduardo Mireles-Cabodevila
Cleveland Clinic Lerner College of Medicine of Case
Western Reserve
Simulation and Advanced Skills Center, Cleveland Clinic
Cleveland, Ohio

Marcelo Gama de Abreu
 Cleveland Clinic Lerner College of Medicine of Case
 Western Reserve
 Department of Intensive Care and Resuscitation,
 Department of Outcomes Research,
 Anesthesiology Institute, Cleveland Clinic
 Cleveland, Ohio

REFERENCES

- Vargas M, Buonanno P, Sica A, Ball L, Iacovazzo C, Marra A, et al. Patient-ventilator synchrony in neurally adjusted ventilator assist (NAVA) and variable pressure support ventilation (PSV). *Respir Care* 2022;67(5):503-509.
- Chatburn RL, El-Khatib M, Mireles-Cabodevila E. A taxonomy for mechanical ventilation: 10 fundamental maxims. *Respir Care* 2014;59(11):1747-1763.
- Mireles-Cabodevila E, Hatipoğlu U, Chatburn RL. A rational framework for selecting modes of ventilation. *Respir Care* 2013;58(2):348-366.
- Mireles-Cabodevila E, Siuba MT, Chatburn RL. A taxonomy for patient-ventilator interactions and a method to read ventilator waveforms. *Respir Care* 2022;67(1):129-148.
- Thille AW, Lyazidi A, Richard J-CM, Galia F, Brochard L. A bench study of intensive care unit ventilators: new versus old and turbine-based versus compressed gas-based ventilators. *Intensive Care Med* 2009;35(8):1368-1376.
- Brochard L, Telias I. Bedside detection of overassistance during pressure support ventilation.* *Crit Care Med* 2018;46(3):488-490.
- Piquilloud L, Vignaux L, Bialais E, Roeseler J, Sottiaux T, Laterre P-F, et al. Neurally adjusted ventilatory assist improves patient-ventilator interaction. *Intensive Care Med* 2011;37(2):263-271.
- Carteaux G, Córdoba-Izquierdo A, Lyazidi A, Heunks L, Thille AW, Brochard L. Comparison between neurally adjusted ventilatory assist and pressure support ventilation levels in terms of respiratory effort. *Crit Care Med* 2016;44(3):503-511.
- Ferreira JC, Diniz-Silva F, Moriya HT, Alencar AM, Amato MBP, Carvalho CRR. Neurally adjusted ventilatory assist (NAVA) or pressure support ventilation (PSV) during spontaneous breathing trials in critically ill patients: a crossover trial. *BMC Pulm Med* 2017;17(1):139.
- Diniz-Silva F, Moriya HT, Alencar AM, Amato MBP, Carvalho CRR, Ferreira JC. Neurally adjusted ventilatory assist vs pressure support to deliver protective mechanical ventilation in patients with acute respiratory distress syndrome: a randomized crossover trial. *Ann Intensive Care* 2020;10(1):18.
- Lamouret O, Crognier L, Vardon Bounes F, Conil J-M, Dilasser C, Raimondi T, et al. Neurally adjusted ventilatory assist (NAVA) versus pressure support ventilation: patient-ventilator interaction during invasive ventilation delivered by tracheostomy. *Crit Care* 2019;23(1):2-9.
- Thille AW, Rodriguez P, Cabello B, Lellouche F, Brochard L. Patient-ventilator asynchrony during assisted mechanical ventilation. *Intensive Care Med* 2006;32(10):1515-1522.
- Rolland-Debord C, Bureau C, Poitou T, Belin L, Clavel M, Perbet S, et al. Prevalence and prognosis impact of patient-ventilator asynchrony in early phase of weaning according to two detection methods. *Anesthesiology* 2017;127(6):989-997.
- Sinderby C, Navalesi P, Beck J, Skrobik Y, Comtois N, Friberg S, et al. Neural control of mechanical ventilation in respiratory failure. *Nat Med* 1999;5(12):1433-1436.
- de la Oliva P, Schüffelmann C, Gómez-Zamora A, Villar J, Kacmarek RM. Asynchrony, neural drive, ventilatory variability, and COMFORT: NAVA versus pressure support in pediatric patients. A nonrandomized crossover trial. *Intensive Care Med* 2012;38(5):838-846.
- Demoule A, Clavel M, Rolland-Debord C, Perbet S, Terzi N, Kouatchet A, et al. Neurally adjusted ventilatory assist as an alternative to pressure support ventilation in adults: a French multi-center randomized trial. *Intensive Care Med* 2016;42(11):1723-1732.
- Haudebourg A-F, Perier F, Tuffet S, de Prost N, Razazi K, Mekontso Dessap A, et al. Respiratory mechanics of COVID-19– versus non-COVID-19–associated acute respiratory distress syndrome. *Am J Respir Crit Care Med* 2020;202(2):287-290.
- Grieco DL, Bitondo MM, Aguirre-Bermeo H, Italiano S, Idone FA, Moccaldò A, et al. Patient-ventilator interaction with conventional and automated management of pressure support during difficult weaning from mechanical ventilation. *J Crit Care* 2018;48:203-210.
- Magrans R, Ferreira F, Sarlabous L, López-Aguilar J, Gomà G, Fernandez-Gonzalo S, et al. The effect of clusters of double triggering and ineffective efforts in critically ill patients. *Crit Care Med* 2022. Publish Ahead of Print.
- Younes M. Why does airway pressure rise sometimes near the end of inflation during pressure support? *Intensive Care Med* 2008;34(1):1-3.
- Doorduyn J, Sinderby CA, Beck J, van der Hoeven JG, Heunks LMA. Assisted ventilation in patients with acute respiratory distress syndrome: lung-distending pressure and patient-ventilator interaction. *Anesthesiology* 2015;123(1):181-190.