



Vios™ Aerosol Delivery System Configurations



	Part #	NDC #	HCPCS #
Standard, LC® Sprint	310F35-LCS	83490-0310-58	E0570, A7005
Standard, LC® Plus	310F83-LC+	83490-0310-67	E0570, A7005
Deluxe, LC® Sprint	310F35-LSD	83490-0310-59	E0570, A7005
Deluxe, LC® Plus	310F81-LCD	83490-0310-64	E0570, A7005



All configurations include Wing Tip™ Tubing





(OR)











Life Inspired

Technical Specifications:

- Maximum Pressure of Compressor = 34 PSI (2.3 bar)
- Maximum Liter Flow of Compressor = 8.5 LPM
- Nominal Operating Pressure² = 18 PSI (1.2 bar)
- Nominal Air Flow² = 4.5 LPM
- Weight = 3 lbs (1300 gms)
- Dimensions (Length, Width, Height) = 6.5 x 6.5 x 3.75 in
- Warranty = 5 years, 6 years user activated

Vios™ Accessories :

- PARI LC Sprint = 023F35
- PARI LC Plus = 022F81
- Replacement Air Filters = 041E4851P12

²When used with a PARI LC® reusable nebulizer.



INTRODUCING



AEROSOL DELIVERY SYSTEM

From the leader in clinically proven aerosol delivery



PARI nebulizers are the most widely used nebulizers by pharmaceutical companies in clinical trials for new medications compared to any other brand.



WHY VIOS™?



“ Because there is a difference in aerosol delivery systems and Vios™ outperforms the others. ”



“ Because with Vios™, I am confident that my patients will receive an effective treatment dose. ”



“ Because my Vios™ treatments are fast, I can spend more time with my grandchildren. ”

The Vios™ Aerosol Delivery System is designed from PARI's long legacy of clinically proven products and consists of two important components: a Vios™ compressor and a PARI nebulizer. Vios™ was inspired by the Latin and Greek words for “Life”. Our goal is to improve quality of “Life” and patient lifestyle through efficient aerosol delivery. The Vios™ delivers consistent particle size and fast treatments resulting in targeted delivery to the lungs.

Consistent Aerosol Delivery

EXCLUSIVE TO PARI, **TIMESTRIP®** IS AN EASY TO USE, VISUAL INDICATOR THAT CLEARLY INSTRUCTS PATIENTS TO REPLACE THEIR PARI REUSABLE NEBULIZERS EVERY 6 MONTHS. NO MORE GUESSWORK OR HASSLE WITH CALENDARS.



- Easy to use and activate
- Great for Patient Reminder programs
- Benefits of replacing PARI nebulizers:
 1. Consistent medication delivery
 2. Fast treatment times
 3. Reduces the risk for bacterial contamination¹

¹Lester, Mary K. Nebulizer Use and Maintenance by Cystic Fibrosis Patients: A Survey Study. Respiratory Care 2004; 49(12): 1504-1508.